

AccuPoint[®]
Advanced



A New Level of Clean

ATP hygiene monitoring that delivers

Simplicity • Audit Compliance • Brand Protection





Is your hygiene program effective?

Good hygiene monitoring forms the core of food safety programs, and can catch problems before they ever reach consumers. Preventing contamination with pathogens or allergens protects public health — and brand integrity.

Up to 35% of food recalls are related to poor sanitation. Many of these may have been prevented if effective monitoring programs were in place¹.

Neogen's AccuPoint® Advanced system (integrated sampler, reader and software) provides an economical, fast and accurate way to monitor cleanliness by measuring adenosine triphosphate (ATP) on surfaces or in liquids.

No matter if your hygiene program is basic, enhanced or advanced, AccuPoint Advanced is ready when you are.

Immediate Action

It is important to know that ATP is not uniformly distributed on surfaces after cleaning. So, it is critical to use a system that

has the ability to both find and recover ATP from surfaces or liquids. Also, important decisions will be made based on those results, so a system that can accurately measure the level of ATP is important.

The AccuPoint Advanced system's scale allows for higher resolution across a broad range. By using Data Manager software to review results, it can help detect slight changes in plant ATP levels.

The AccuPoint Advanced system integrates the best in sampler design and chemistry, fast time to result and the ability to establish baselines over time to help make informed decisions.

Design Simplicity

AccuPoint Advanced samplers were designed to use minimal plastic. Also, part of the design was to ensure that all recovered ATP gets converted to signal, so having the needed reagents in close proximity helps eliminate any loss of signal.



1. <http://www.webbertraining.com/files/library/docs/32.pdf>.



AccuPoint Advanced Samplers

Since ATP is not evenly distributed on surfaces, Neogen's Advanced sampling pad and chemistry are designed to recover and detect small amounts of ATP.




- Neogen's proprietary sampler is designed specifically to maximize recovery of ATP from surfaces and liquids.
- Newly developed liquid stable chemistry provides more consistent results.
- Samplers can be kept out of refrigeration for up to 2 weeks.
- Internal studies show up to 10 times greater recovery than traditional ATP swabs.

AccuPoint Advanced Instrument

AccuPoint Advanced is the most compact instrument on the market. With updated firmware and RFID technology, it offers the latest technology in the smallest instrument.

- Faster time to result (<20 seconds) allows more data collection in less time.
- RFID technology (optional feature) streamlines the process of collecting hygiene data by choosing the site to be tested based on your preference.
- Unique relative light unit (RLU) scale offers higher resolution across all levels to identify deviations in ATP baselines.



AccuPoint Advanced Data Manager

AccuPoint Advanced's Data Manager software allows users to streamline the testing process by quickly creating test plans, syncing important data, and keeping a permanent record of hygiene results to verify SSOPs.

- View pass rate percentages for entire facility, groups, sites or users to help drive continuous improvement.
- Use the rank report to quickly identify areas that need attention.
- Use the trend report to view RLU baselines.
- Maintain records for audit compliance.
- Firmware upgrade feature keeps reader current.



Corporate Profile

Founded in 1982, Neogen Corporation (NASDAQ: NEOG) has grown to have a worldwide presence, with more than 1,000 employees in multiple U.S. and international locations. Neogen remains firmly committed to consistently meeting the quality standards customers expect.

Neogen's Food Safety Division develops and markets solutions to simply, rapidly and accurately detect bacterial pathogens, mycotoxins, spoilage microorganisms, food allergens, genetic modifications, drug and pesticide residues and hygiene concerns in food and animal feed.



Step up to the complete ATP monitoring system

Ask about a FREE 30 day evaluation.



Informed Decisions

AccuPoint Advanced provides the information users need when they need it.

Prevention and Compliance

AccuPoint Advanced can help producers verify their hygiene and preventive control programs are effective.

Brand Protection

AccuPoint Advanced provides confidence in results that maintain brand reputation.



800-234-5333 (USA/Canada) • 517-372-9200
foodsafety@neogen.com • www.neogen.com

