



W-L NUTRIENT MEDIUM (7488)

Intended Use

W-L Nutrient Medium is used for the cultivation of yeasts, molds, and bacteria encountered in brewing and industrial fermentations in a laboratory setting. W-L Nutrient Medium is not intended for use in the diagnosis of disease or other conditions in humans.

Product Summary and Explanation

W-L Nutrient Medium was developed by Green and Gray^{1,2} while studying various fermentation processes. An exhaustive study examining methods of fermentation control procedures in worts, beers, liquid yeasts and similar fermentation products led to the development of W-L Nutrient Medium. At a pH of 5.5, counts of viable baker's yeast will grow on W-L Nutrient Medium.

W-L Nutrient Medium is also referred to as "Wallerstein Laboratory Medium".

Principles of the Procedure

Yeast Extract is a source of trace elements, vitamins, and amino acids. Enzymatic Digest of Casein provides nitrogen, amino acids, and carbon. Dextrose is a source of carbohydrate. Monopotassium Phosphate buffers the medium. Potassium Chloride, Calcium Chloride, and Ferric Chloride are essential ions and help to maintain osmotic balance. Magnesium Sulfate and Manganese Sulfate are sources of divalent cations. Bromcresol Green is a pH indicator. Agar is the solidifying agent.

Formula / Liter

Yeast Extract.....	4 g
Enzymatic Digest of Casein	5 g
Dextrose.....	50 g
Monopotassium Phosphate	0.55 g
Potassium Chloride.....	0.425 g
Calcium Chloride.....	0.125 g
Magnesium Sulfate	0.125 g
Ferric Chloride	0.0025 g
Manganese Sulfate	0.0025 g
Bromcresol Green.....	0.022 g
Agar	20 g

Final pH: 5.5 ± 0.2 at 25°C

Formula may be adjusted and/or supplemented as required to meet performance specifications.

Precaution

1. For Laboratory Use Only.

Directions

1. Dissolve 80 g of the medium in one liter of purified water.
2. Heat with frequent agitation to completely dissolve the medium.
3. Autoclave at 121°C for 15 minutes.

Quality Control Specifications

Dehydrated Appearance: Powder is homogeneous, free flowing, and light blue or green beige.

Prepared Appearance: Prepared medium is trace to slightly hazy and blue to blue-green.

Expected Cultural Response: Cultural response on W-L Nutrient Medium incubated aerobically at 35 ± 2°C and examined for growth after 18 - 48 hours.

Microorganism	Approx. Inoculum (CFU)	Response
<i>Escherichia coli</i> ATCC® 25922	10 - 300	Growth
<i>Lactobacillus fermentum</i> ATCC® 9338	10 - 300	Growth
<i>Proteus mirabilis</i> ATCC® 12453	10 - 300	Growth
<i>Saccharomyces cerevisiae</i> ATCC® 9763	10 - 300	Growth

The organisms listed are the minimum that should be used for quality control testing.

Test Procedure

Refer to appropriate references for specific procedures. For a complete discussion on the isolation and identification of yeasts, refer to references outlined in the references.^{3,4}

Results

Refer to appropriate references and procedures for results.

Storage

Store sealed bottle containing the dehydrated medium at 2 - 30°C. Once opened and recapped, place container in a low humidity environment at the same storage temperature. Protect from moisture and light by keeping container tightly closed.

Expiration

Refer to expiration date stamped on the container. The dehydrated medium should be discarded if not free flowing, or if appearance has changed from the original color. Expiry applies to medium in its intact container when stored as directed.

Limitation of the Procedure

Due to varying nutritional requirements, some strains may be encountered that grow poorly or fail to grow on this medium.

Packaging

W-L Nutrient Medium	Code No.	7488A	500 g
		7488B	2 kg
		7488C	10 kg

References

1. **Green, S. R., and P. P. Gray.** 1950. Paper read at American Society of Brewing Chemists Meeting. Wallerstein Lab. Commun. 12:43.
2. **Green, S. R., and P. P. Gray.** 1950. A differential procedure applicable to bacteriological investigation in brewing. Wallerstein Lab. Commun. 13:357.
3. **Murray, P. R., E. J. Baron, M. A. Pfaller, F. C. Tenover, and R. H. Tenover (eds.).** Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D. C.
4. **Isenberg, H. D. (ed.).** 1992. Interpretation of aerobic bacterial growth on primary culture media, Clinical microbiology procedures handbook, vol. 1 p. 1.61-1.6.7. American Society for Microbiology, Washington, D.C.

Technical Information

Contact Acumedia Manufacturers, Inc. for Technical Service or questions involving dehydrated culture media preparation or performance at (517)372-9200 or fax us at (517)372-2006.