



Agri-Screen Ticket




Intended Use

Agri-Screen Ticket is used to detect all major organophosphates, thiophosphates and carbamates in air, water, soil, produce, food, spills, surfaces and custom applications. Among the common pesticides that Ticket detects are: Malathion, Sevin, Furadan, Systox, Vapona, Phosdrin, Metasystox, Guthion, Actellic, Dursban and diazinon. For a complete list, contact Neogen. Samples can be tested wet, dry, or with solvents.

Results



Positive



Negative

Ordering Information

Prod.#	Product description
8920	Agri-Screen Ticket-10 Kit – Contains 10 tickets, 10 ampoules, beaker and rod
8920A	Agri-Screen Ticket-10 Kit – Contains 10 tickets, beaker and rod
8921	Agri-Screen Ticket-25 Kit – Contains 25 tickets, 25 ampoules, beaker and rod
8921A	Agri-Screen Ticket-25 Kit – Contains 25 tickets, beaker and rod
8950	Agri-Screen Tickets – Contains 10 tickets
8951	Activator Ampoules – Contains 10 ampoules
8956	Activator Ampoules – Contains 5 ampoules
8953	Concentrator
8954	Beaker
8955	Glass rod

Neogen Corporation

800/234-5333 • 517/372-9200 • E-mail: foodsafety@neogen.com

Neogen Europe

+44 (0) 1292 525600 • E-mail: info@neogeneurope.com