



Neogen 4700 Reader (9303)

- Automatic strip transportation
- Bichromic optics with four filters:
Neogen (Standard) 405, 450, 492 and 650 nm
- Easy to use with step-by-step prompting
- On-board graphics printing
- Calculates results automatically, large nonvolatile memory stores more than 100 user-entered tests
- Includes one 3 x 12 strip carrier for break apart microwells
- Syncs with Veratox software for simple reporting, storage and tracking of results

Neogen Corporation

800/234-5333 • 517/372-9200 • E-mail: foodsafety@neogen.com

Neogen Europe

+44 (0) 1292 525275 • E-mail: info@neogeneurope.com