

# Reveal

## For CP4 (Roundup Ready)




### Intended Use

Reveal for CP4 is an immunochromatographic lateral flow assay used for the qualitative analysis of CP4 in soybeans and corn. The test can detect as little as 1 CP4 enhanced soybean among 1,000 beans (0.1%), and 1 NK603 Roundup Ready kernel of corn among 800 kernels (0.125%). Earlier Roundup Ready corn varieties, such as GA21, will not be detected.

### Product Specifications

Sensitivity:	Soybeans: 1 in 1,000 beans (0.1%) Corn: 1 in 800 kernels (0.125%)
Testing time:	5–10 minutes
Storage:	Room temperature (18–30°C)
Tests per kit:	Up to 50

### Results



### Ordering Information

Prod.#	Product description
8005	Reveal for CP4 – screens 50 samples

### Neogen Corporation

800/234-5333 • 517/372-9200 • E-mail: [foodsafety@neogen.com](mailto:foodsafety@neogen.com)

### Neogen Europe

+44 (0) 1292 525600 • E-mail: [info@neogeneurope.com](mailto:info@neogeneurope.com)