

Reveal 2.0

For *Salmonella*




Intended Use

The Reveal 2.0 for *Salmonella* test provides for the rapid recovery of *Salmonella* in chicken carcass rinse, raw ground turkey, raw ground beef, hot dogs, raw shrimp, ready-to-eat meal products, dry pet food, ice cream, fresh spinach, cantaloupe, peanut butter, swabs from stainless steel surfaces, and sprout irrigation water samples allowing detection and presumptive identification of *Salmonella* within 24 hours.

Product Specifications

Sensitivity:	1 cfu/analytical unit 10 ⁶ cfu/mL post enrichment
Testing time:	15 minutes
Tests per kit:	20
Kit contains:	20 Reveal 2.0 for <i>Salmonella</i> test devices 20 disposable transfer pipettes 20 sample cups
Approvals:	AOAC-RI 960801

Results



Ordering Information

Prod.#	Product description
9706	Reveal 2.0 One-Step for <i>Salmonella</i> – Kit contains 20 devices
9802	Reveal 2.0 One-Step <i>Salmonella</i> Complete System – Contains devices, RV selective medium, and sampling bags for 20 tests
9803	Reveal 2.0 for <i>Salmonella</i> Complete System with RV – Contains devices, REVIVE, RV selective medium and sampling bags for 20 tests
9804	Reveal 2.0 for <i>Salmonella</i> Complete System with M-Broth and RV – Contains devices, REVIVE, RV selective medium, M-Broth and sampling bags for 20 tests
9805	Reveal 2.0 for <i>Salmonella</i> Environmental Sampling System – Contains devices, REVIVE, RV selective medium, M-Broth and environmental sampling kits for 20 tests
9705	REVIVE for <i>Salmonella</i> – Patented resuscitation medium to speed up enrichment for <i>Salmonella</i> . Contains bottles and bags for 20 tests
9729	1X RV Selective Media – Rappaport-Vassiliadis Broth 1X. Media kit contains prefilled bottles for 20 tests and 20 stomacher bags
9715	2X RV Selective Media – Rappaport-Vassiliadis Broth 2X concentrate. Media kit contains prefilled bottles for 20 tests and 20 stomacher bags. To be used in combination with REVIVE.

Neogen Corporation

800/234-5333 • 517/372-9200 • E-mail: foodsafety@neogen.com

Neogen Europe

+44 (0) 1292 525600 • E-mail: info@neogeneurope.com