

Soleris

For Generic *E. coli* and Total Coliforms



Intended Use

Soleris *E. coli* and Coliform Vials are used as a screening method for detection of a variety of coliform organisms.

Product Specifications



Testing time:	4–24 hours
Tests per kit:	100
Approvals:	AOAC-RI 101101

Materials Required But Not Provided

Available from Neogen

1. Soleris Automated System
2. Sterile PBS dilution bottles (BPB-99)
3. Adjustable pipette 1 mL
4. Tryptone (Indole) broth confirmation (TT-123)
5. Kovac's reagent (KO-122)
6. Tryptic Soy Broth (TS-124)

Results



Ordering Information

Prod.#	Product description
CC-102	Coliform Medium 5 mL (AOAC-RI)
CC-109	Coliform Medium 9 mL (AOAC-RI)
CC-BP9	Coliform Medium Buffered 9 mL
CC-BP9-PRO	Coliform Medium Buffered 9 mL, low pH
EC-104	<i>E. coli</i> Medium (AOAC-RI)

Neogen Corporation

800/234-5333 • 517/372-9200 • E-mail: foodsafety@neogen.com

Neogen Europe

+44 (0) 1292 525600 • E-mail: info@neogeneurope.com