

Veratox

For Aflatoxin M₁




Intended Use

Veratox for Aflatoxin M₁ is intended for the quantitative analysis of aflatoxin M₁ in commodities such as milk, milk powder, butter and cheese.

Product Specifications

Lower limit of detection:	4.3 ppt
Range of quantitation:	5–100 ppt
Controls provided:	0, 5, 15, 30, 60 and 100 ppt
Testing time:	45 minutes
Antibody cross-reactivity:	Aflatoxin M ₁ 100%. Aflatoxin B ₁ , B ₂ , G ₁ , G ₂ and M ₂ < 1%.
Tests per kit:	Up to 38

Results



Ordering Information

Prod.#	Product description
8019	Veratox for Aflatoxin M ₁

Neogen Corporation

800/234-5333 • 517/372-9200 • E-mail: foodsafety@neogen.com

Neogen Europe

+44 (0) 1292 525600 • E-mail: info@neogeneurope.com