

Veratox

For Crustacea Allergen




Intended Use

Veratox for Crustacea Allergen is used for the quantitative analysis of crustacea proteins in food products, ingredients and clean-in-place rinses.

Product Specifications

Range of quantitation:	2.5–25 ppm
Controls provided:	0, 2.5, 5, 10 and 25 ppm crustacea
Testing time:	30 minutes
Tests per kit:	Up to 38

Results



Ordering Information

Prod.#	Product description
8520	Veratox for Crustacea Allergen
8429	Allergen Extraction Kit – 20 samples
8432S	Allergen Environmental Swabbing Kit – 100 swabs and tips

Neogen Corporation

800/234-5333 • 517/372-9200 • E-mail: foodsafety@neogen.com

Neogen Europe

+44 (0) 1292 525600 • E-mail: info@neogeneurope.com