

# Veratox®

## For Sesame Allergen




### Intended Use

Veratox for Sesame Allergen is used for the quantitative analysis of sesame residue in food products such as spices, bakery products, cereal products and clean-in-place rinses.

### Product Specifications

Range of quantitation:	2.5–25 ppm
Controls provided:	0, 2.5, 5, 10 and 25 ppm sesame
Testing time:	30 minutes
Tests per kit:	Up to 38

### Results



### Ordering Information

Prod#.	Product description
8530	Veratox for Sesame Allergen
8429	Allergen Extraction Kit - 20 samples

### Neogen Corporation

800/234-5333 • 517/372-9200 • E-mail: [foodsafety@neogen.com](mailto:foodsafety@neogen.com)

### Neogen Europe

+44 (0) 1292 525600 • E-mail: [info@neogeneurope.com](mailto:info@neogeneurope.com)