

B-Pep CCT/ET 01 (NCM0261)

Intended Use

A blended tryptose derived from enzymatic digest of animal origin meat and casein peptones developed as a nutritional supplement for use in animal cell culture applications. B-Pep CCT/ET 01 is not intended for use in the diagnosis of disease or other conditions in humans.

Description

A blended tryptose derived from enzymatic digest of animal origin meat and casein peptones developed as a nutritional supplement for use in animal cell culture applications.

Precaution

Refer to SDS

Test Procedure

Refer to appropriate references for specific procedures using B-Pep CCT/ET 01.

Quality Control Specifications

PHYSICAL/CHEMISTRY:	EXPECTED RESULTS
Appearance of Powder:	Fine, free-flowing, light tan powder
Clarity:	Clear, yellow color
Solubility In water:	Total
pH of 2% solution	7.0 – 7.4
Endotoxin Assay (20 g/L)	<4 EU/mL

Expiration

Refer to expiration date stamped on container. B-Pep CCT/ET 01 should be discarded if not free flowing, or if the appearance has changed from the original color. Expiry applies to B-Pep CCT/ET 01 in its intact container when stored as directed.

Storage

Store powder at 2-30°C away from direct sunlight. Once opened and recapped, place container in a low humidity environment at the same storage temperature. Protect from moisture and light by keeping container tightly closed.