

## M-Pep CCT/ET 01 (NCM0262)

### Intended Use

An enzymatic digest of meat, this peptone is highly nutritious and has been developed as a nutritional supplement for use in animal cell culture applications. M-Pep CCT/ET 01 is not intended for use in the diagnosis of disease or other conditions in humans.

### Description

An enzymatic digest of meat, this peptone is highly nutritious and has been developed as a nutritional supplement for use in animal cell culture applications

### Precaution

Refer to SDS

### Test Procedure

Refer to appropriate references for specific procedures using M-Pep CCT/ET 01.

### Quality Control Specifications

<b>PHYSICAL/CHEMISTRY:</b>	<b>EXPECTED RESULTS</b>
Appearance of Powder:	Fine, free-flowing, beige powder
Clarity:	Clear, yellow color
Solubility In water:	Total
pH of 2% solution	6.5 – 7.5
Endotoxin Assay (20 g/L)	<4 EU/mL

### Expiration

Refer to expiration date stamped on container. M-Pep CCT/ET 01 should be discarded if not free flowing, or if the appearance has changed from the original color. Expiry applies to M-Pep CCT/ET 01 in its intact container when stored as directed.

### Storage

Store powder at 2-30°C away from direct sunlight. Once opened and recapped, place container in a low humidity environment at the same storage temperature. Protect from moisture and light by keeping container tightly closed.