

Agri-Screen

For Fumonisin





Intended Use

Agri-Screen for Fumonisin is a competitive direct ELISA that provides a visible screening result against a known control of fumonisin in such commodities as barley, corn, DDGS, milo, popcorn, rice, soybeans and wheat.

Product Specifications

Control provided:	5 ppm
Testing time:	15 minutes
Antibody cross-reactivity:	Total fumonisins (B ₁ , B ₂ , B ₃)
Test per kit:	Up to 20

Results



Ordering Information

Prod.#	Product description
8810	Agri-Screen for Fumonisin
8873	Fumonisin 5/10 Dilution Bottles – 28 samples
8052	Mycotoxin Extraction Kit

Neogen Corporation

800/234-5333 • 517/372-9200 • E-mail: foodsafety@neogen.com

Neogen Europe

+44 (0) 1292 525600 • E-mail: info@neogeneurope.com