


# Biokits F.A.S.T. Immunostick Meat Species Screening Kit




## Intended Use

Biokits F.A.S.T. Immunostick Meat Species Screening Kit is intended for the qualitative detection of species content in raw, processed and mixed meats and milk.

## Product Specifications

Sensitivity:	Screens at 1%
Testing time:	30 minutes
Storage:	Refrigerated 2–8°C
Tests per kit:	25 (5 sets of 5 test packs)

## Results



## Ordering Information

Prod.#	Product description
502135B	F.A.S.T. Cow – Qualitative Immunostick tube assays capable of detecting 1% of animal species. (5 tests per pack)
502136Y	F.A.S.T. Horse – Qualitative Immunostick tube assays capable of detecting 1% of animal species. (5 tests per pack)
502137W	F.A.S.T. Pig – Qualitative Immunostick tube assays capable of detecting 1% of animal species. (5 tests per pack)
502139S	F.A.S.T. Poultry – Qualitative Immunostick tube assays capable of detecting 1% of animal species. (5 tests per pack)
502138U	F.A.S.T. Sheep – Qualitative Immunostick tube assays capable of detecting 1% of animal species. (5 tests per pack)
502304E	F.A.S.T. Positive Controls

## Neogen Corporation

800/234-5333 • 517/372-9200 • E-mail: [foodsafety@neogen.com](mailto:foodsafety@neogen.com)

## Neogen Europe

+44 (0) 1292 525600 • E-mail: [info@neogeneurope.com](mailto:info@neogeneurope.com)