

NeoColumn

For Ochratoxin A (HPLC Clean-up)



Intended Use

NeoColumn for Ochratoxin A is a high performance immunoaffinity column designed for clean-up and concentration of sample prior to running on an HPLC system or Veratox. NeoColumn is suitable for use in the determination of ochratoxin A in a wide range of commodities such as cereals, dried fruits, beer and wine, and green and roast coffee.


Product Specifications

Lower limit of detection:	< 0.1 ppb
Testing time:	30 minutes
Recovery:	> 95%*
Columns per kit:	50

*Example given is for cereals. Conditions may vary depending on the commodity.

Results

Chromatographic detection
of ochratoxin A



Ordering Information

Prod.#	Product description
8640	NeoColumn for Ochratoxin A
8651	Extraction Kit for Wine and Grape Juice – 25 samples
8652	Wine/Grape Juice Clarification Compound – 50 packs
9371	Syringe adapter for NeoColumn
9379	5 mL syringe for NeoColumn

Neogen Corporation

800/234-5333 • 517/372-9200 • E-mail: foodsafety@neogen.com

Neogen Europe

+44 (0) 1292 525600 • E-mail: info@neogeneurope.com