

# Soleris

## Direct Vial for Yeast and Mold



### Intended Use

Soleris Yeast and Mold Direct Vials are used as a screening method for detection of a variety of yeast and mold spoilage.

### Product Specifications

Testing time: 4–48 hours


Tests per kit: 100

### Materials Required But Not Provided

Available from Neogen

1. Soleris Automated System
2. Sterile PBS dilution bottles
3. Yeast and mold supplement (YI-110)

### Results



### Ordering Information

Prod.#	Product description
DYM-109	Direct Vial for Yeast and Mold
YI-110	Yeast and Mold Supplement

### Neogen Corporation

800/234-5333 • 517/372-9200 • E-mail: [foodsafety@neogen.com](mailto:foodsafety@neogen.com)

### Neogen Europe

+44 (0) 1292 525600 • E-mail: [info@neogeneurope.com](mailto:info@neogeneurope.com)