

# Veratox

## For Fumonisin 5/10




### Intended Use

Veratox for Fumonisin 5/10 is a competitive direct ELISA that provides a quantitative analysis of fumonisin in commodities such as barley, corn, DDGS, milo, popcorn, rice, soybeans and wheat.

### Product Specifications

Lower limit of detection:	0.2 ppm
Range of quantitation:	0.5–6 ppm
Controls provided:	0, 0.5, 1, 3 and 6 ppm
Testing time:	15 minutes
Antibody cross-reactivity:	Total fumonisins (B <sub>1</sub> , B <sub>2</sub> , B <sub>3</sub> )
Tests per kit:	Up to 38

### Results



### Ordering Information

Prod.#	Product description
8835	Veratox for Fumonisin 5/10
8873	Fumonisin 5/10 Dilution Bottles – 28 samples
8052	Mycotoxin Extraction Kit

### Neogen Corporation

800/234-5333 • 517/372-9200 • E-mail: [foodsafety@neogen.com](mailto:foodsafety@neogen.com)

### Neogen Europe

+44 (0) 1292 525600 • E-mail: [info@neogeneurope.com](mailto:info@neogeneurope.com)