

# Veratox

## For Fumonisin




### Intended Use

Veratox for Fumonisin is a competitive direct ELISA that provides a quantitative analysis of fumonisin in commodities such as barley, corn, DDGS, milo, popcorn, rice, soybeans and wheat.

### Product Specifications

Lower limit of detection:	0.2 ppm
Range of quantitation:	1–6 ppm
Controls provided:	0, 1, 2, 4 and 6 ppm
Testing time:	20 minutes
Antibody cross-reactivity:	Total fumonisins (B <sub>1</sub> , B <sub>2</sub> , B <sub>3</sub> )
Tests per kit:	Up to 38
Approvals:	AOAC Official Method 2001.06

### Results



### Ordering Information

Prod.#	Product description
8830	Veratox for Fumonisin with dilution
8872	Fumonisin Dilution Bottles – 14 samples
8052	Mycotoxin Extraction Kit

### Neogen Corporation

800/234-5333 • 517/372-9200 • E-mail: [foodsafety@neogen.com](mailto:foodsafety@neogen.com)

### Neogen Europe

+44 (0) 1292 525600 • E-mail: [info@neogeneurope.com](mailto:info@neogeneurope.com)