

Veratox

For Fumonisin HS (high sensitivity)




Intended Use

Veratox for Fumonisin HS is a competitive direct ELISA that provides a quantitative analysis of fumonisin in commodities such as barley, corn, DDGS, milo, popcorn, rice, soybeans and wheat.

Product Specifications

Lower limit of detection:	50 ppb
Range of quantitation:	50–600 ppb
Controls provided:	0, 50, 100, 300 and 600 ppb
Testing time:	15 minutes
Antibody cross-reactivity:	Total fumonisins (B ₁ , B ₂ , B ₃)
Tests per kit:	Up to 38

Results



Ordering Information

Prod.#	Product description
8832	Veratox for Fumonisin HS
8052	Mycotoxin Extraction Kit

Neogen Corporation

800/234-5333 • 517/372-9200 • E-mail: foodsafety@neogen.com

Neogen Europe

+44 (0) 1292 525600 • E-mail: info@neogeneurope.com