

Veratox®

For Lincomycin




Intended Use

Veratox for Lincomycin is used for the quantitative analysis of lincomycin in muscle.

Product Specifications

Assay Range: 7 - 262.5 ppb
Controls provided: 0, 0.1, 0.25, 0.6, 1.5 and 3.75 ng/mL Lincomycin
Testing time: 50 minutes
Tests per kit: Up to 36 run in duplicate

Quantitative Results



Ordering Information

Prod.#	Product description
DR206	Veratox for Lincomycin