

# Veratox HS

## For Ochratoxin




### Intended Use

Veratox HS for Ochratoxin is a competitive direct ELISA that provides quantitative analysis of lower levels of ochratoxin in wheat and corn.

### Product Specifications

|                           |                       |
|---------------------------|-----------------------|
| Lower limit of detection: | 1 ppb                 |
| Range of quantitation:    | 2–10 ppb              |
| Controls provided:        | 0, 3, 7.5, and 15 ppb |
| Testing time:             | 30 minutes            |
| Tests per kit:            | Up to 38              |

### Results



### Ordering Information

| Prod.# | Product description       |
|--------|---------------------------|
| 8632   | Veratox HS for Ochratoxin |
| 8052   | Mycotoxin Extraction Kit  |

### Neogen Corporation

800/234-5333 • 517/372-9200 • E-mail: [foodsafety@neogen.com](mailto:foodsafety@neogen.com)

### Neogen Europe

+44 (0) 1292 525600 • E-mail: [info@neogeneurope.com](mailto:info@neogeneurope.com)